

---



# PEARLCATCH™ SINGLE CHAMBER POLYP TRAP

---

*Instructions for Use.  
Read carefully prior to use.*



**GA Health Company Limited**  
Unit 18, 21/F, Metropole Square,  
2 On Yiu Street, Shatin, N.T., Hong Kong  
Tel +852 2833 9010  
www.GAhealth.net  
info@GAHealth.net

***Not made with natural rubber latex.  
Does not contain DEHP***

<u>Part Number</u>	<u>Description</u>
GAR014	PEARLCATCH™ Single Chamber Polyp Trap (Non-sterile)

## **INTENDED USE**

The single use **PearlCatch™ Single Chamber Polyp Trap** is used for suction retrieval and transportation of endoscopically removed polyps for histological examination. The Single Chamber Polyp Trap has one chamber to collect the polyps and one open area for straight suction.

The single use Single Chamber Polyp Trap consists of:

- One (1) jar with one (1) filtered trap to collect polyps, and one (1) open areas for straight suction.
- A lid with one tube for connection to the endoscope.
- One Polyp Pick.

## **WARNINGS AND PRECAUTIONS**

Products marked “single-use” are for single-use only. Do not reuse, reprocess, or resterilize single-use products. The materials used in the manufacture of the device may not withstand repeated reprocessing. The device may not perform as intended by the manufacturer if it is reused. This may lead to failure of the device to perform as intended and/or material degradation, which may result in patient injury, illness or death. Reuse, reprocessing, or resterilization may also increase the risk of contamination of the device and/or cause patient infection or cross-infection, including but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Do not disconnect the pre-installed tubing from the polyp trap lid outer port as this may result in damage to the polyp trap lid and may result in exposure to the contents of the trap during removal.

Protective garments (i.e. gloves, protective eyewear, gowns, etc) should be worn during the use of this device to protect against the risk of cross- contamination and injury.

## **PRODUCT DISPOSAL**

After use, this product may be a potential biohazard. Handle and dispose of the product in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

## **INSTALLING THE POLYP TRAP**

### **Connection:**

1. Connect barbed fitting of trap to suction system.
2. Connect tube side of trap to endoscope suction fitting.

### **Capture:**

3. Apply suction via the endoscope suction valve to trap polyp.

**Polyp Removal and Transport:**


4. After the polyp is trapped, remove the trap from the endoscope and suction system.
5. Separate trap by rotating both halves in opposite directions.
6. Remove polyp and place in specimen jar with fixative, or place trap in specimen jar and fill with fixative. **Note:** *To remove the retrieved polyp, use the provided Polyp Pick to transfer the polyp from the strainer to the appropriate specimen container for histopathology in accordance with institutional guidelines.*
7. Transport to pathology.

**WARNING**
















An issued or revision date for these instructions is included for the user's information. In the event three years has elapsed between this date and product use, the user should contact GA Health Company Limited to determine if additional information is available. Report any serious incident that has occurred in relation to the device to the distributor or manufacturer or competent authority of the member state.

ANDORATE® is a registered trademark of GA Health Company Limited.

PEARLCATCH™ is a registered trademark of GA Health Company Limited.



**EXPLANATION OF SYMBOLS USED ON LABELS AND INSTRUCTIONS FOR USE**

	Use-by date
	Batch or lot code
	Manufacturer
	Date of Manufacture
	Medical Device
	Catalogue number Product Code
	EC Conformity Quality Assurance
	Do not re-use Single use
	Product not made with natural rubber latex
	Do not contain DEHP
	Consult Instructions for Use
	Authorized representative in the European Community
	Keep dry
	Keep away from sunlight
	Federal law (U.S.A.) restricts this device to sale, distribution and use by or on the order of a physician.