



Valve Storage Bag w/ Drawstring - Directions for Use

Order # ZUTR50528

Description

This tough, lightweight polyester mesh bag helps ensure the traceability of critical endoscope components after reprocessing, by keeping the valves and endoscope together. The valve storage bag with drawstring is designed to keep valves and accessories that have been removed from endoscope during manual cleaning, to ensure decontamination of all surfaces, together with the endoscope. The bag is 3" by 5" and features an easy to use drawstring closure with a barrel style cord lock. Compatible with high level disinfectants, gas plasma and EtO sterilization.

Instructions

1. Place the endoscope valves and accessories inside the mesh bag.
2. Holding the cord, gather the top of the bag to close the top.
3. Press the button on the toggle and pull it along to the closed bag and hold it in place.
4. Secure bag to endoscope.
5. Bag, valves, and accessories are to remain with endoscope for tracking purposes.

Specifications

- Order # **ZUTR50528**
- Qty 10
- Tough, lightweight
- White mesh drawstring bag