


Zutron Medical is proud to announce the **Endo-Flush** endoscope flushing aid



Product Specifications

Dimensions:

4.50"(L) x 7.20"(W) x 8.30"(H)

Weight:

10 pounds

Timer Control:

Digital Display with touch keypad

Flush Time:

Default: 2 minutes

Min/Max: 15 seconds - 3 minutes

Electrical:

Input 24 VDC 3.3 amp positive center pole

Warranty:

One year

Endo-Flush Features

- Delivers detergent quickly**
 Provides fluid delivery throughout every endoscope channel quickly and effectively.
- Eases the manual cleaning process**
 Hands free operation reduces fatigue, strain, and repetitive action injuries of your employees versus manually syringe flushing.
- Endo-Flush Display**
 Facility customizable adjustment of flush time and other settings. Visual output reading for ease of adjustment to meet scope time frames. Audible alert once the flush cycle is complete.
- Compact unit**
 Unit may be placed on countertop or wall mounted with the included bracket.

Endo-Flush Ordering Information

Order #	Description
ZUTR30004	Endo-Flush

*Endo-Flush includes wall mount brackets and tubing set

Additional Ordering Information

Order #	Description
ZUTR30001	Complete Tubing Set
ZUTR30001-1	Pick-Up Tube
ZUTR30001-2	2-Way Flush Tube
ZUTR30001-3	3-Way Flush Tube
ZUTR30001-4	4-Way Flush Tube
ZUTR30001-5	Flush Tube Extension
ZUTR30001-6	Flow Verification Meter
ZUTR30001-7	Pick-Up Tube Extension

www.zutron.com

ZM[®]
ZUTRON MEDICAL

877-343-5873