


Zutron Medical is proud to announce the **Endo-Flush** endoscope flushing aid



Product Specifications

Dimensions:

4.50"(L) x 7.20"(W) x 8.30"(H) Case
5.60"(L) x 7.30"(W) x 10.50"(H) Maximum

Weight:

10 pounds

Timer Control:

Digital Display with touch keypad

Flush Time Range:

15 seconds - 3 minutes

Electrical:

Input 24 VDC 3.3 amp positive center pole

Warranty:

One year

Endo-Flush Features

- Delivers detergent quickly**
Provides fluid delivery throughout every endoscope channel quickly and effectively.
- Eases the manual cleaning process**
Hands free operation reduces fatigue from manually syringing scopes and reduces strain.
- Endo-Flush Display**
Allows user to adjust flush times and other settings. Audible alert once the flush cycle is complete.
- Compact unit**
Unit may be placed on countertop or wall mounted with supplied brackets.

Endo-Flush Ordering Information

Order #	Description
ZUTR30004	Endo-Flush

*Endo-Flush includes wall mount brackets and tubing set

Accessory Ordering Information

Order #	Description
ZUTR30001-2	Additional 2-Way Flush Tube
ZUTR30001-3	Additional 3-Way Flush Tube
ZUTR30001-4	Additional 4-Way Flush Tube

Cleanser Ordering Information

Order #	Description
ZUTR28601	Non-Foaming Multi-Enzymatic

www.zutron.com



877-343-5873