

# Zutron Medical would like to make your day easier with our Hand Held Endoscope Leak Testers



## Product Features

- Effective in finding endoscope leaks**  
 Utilization of the incorporated gauge allows for visualizing a drop in pressure to determine if there is a leak.
- Fast and easy**  
 Simple to use and quickly pressurizes the endoscope.
- Compatible with all flexible scopes**  
 Available for FUJIFILM, ACMI, Olympus, Pentax, Richard Wolf, and Storz endoscopes.
- Quickly deflates endoscope after use**  
 Immediately deflate the endoscope with the quick release button once the leak test is finished.

## Hand Held Leak Tester Ordering Information

Order #	Description
ZUTR10002F	Hand Held Endoscope Leak Tester - FUJIFILM
ZUTR10002A	Hand Held Endoscope Leak Tester - ACMI
ZUTR10002O	Hand Held Endoscope Leak Tester - Olympus
ZUTR10002P	Hand Held Endoscope Leak Tester - Pentax
ZUTR10002RW	Hand Held Endoscope Leak Tester - Richard Wolf
ZUTR10002S	Hand Held Endoscope Leak Tester - Storz

## Additional Leak Testers

Order #	Description
ZUTR30005KIT	Endo-Leak Tester Kit
ZUTR30005	Endo-Leak Tester
ZUTR10003	Automatic Leak Tester