

Zutron Medical would like to make your day easier with our **Hand Held Endoscope Leak Testers**



Product Features

- Effective in finding endoscope leaks**
 Watch the needle for a drop in pressure to determine if there is a leak in the endoscope.
- Fast and easy**
 Simple to use and quickly pressurizes the endoscope.
- Compatible with all flexible scopes**
 The Hand Held Endoscope Leak Testers are available for ACMI, Fujinon, Olympus, Pentax, Richard Wolf, and Storz endoscopes.
- Quickly deflates endoscope after use**
 The quick release button immediately deflates the endoscope once the leak test is finished.

Hand Held Leak Tester Ordering Information

Order #	Description
ZUTR10002A	Hand Held Endoscope Leak Tester - ACMI
ZUTR10002F	Hand Held Endoscope Leak Tester - Fujinon
ZUTR10002O	Hand Held Endoscope Leak Tester - Olympus
ZUTR10002P	Hand Held Endoscope Leak Tester - Pentax
ZUTR10002RW	Hand Held Endoscope Leak Tester - Richard Wolf
ZUTR10002S	Hand Held Endoscope Leak Tester - Storz

Accessories Ordering Information

Order #	Description
ZUTRAW	Olympus ETO Alignment Wrench
SEHWM	Single Scope Wall Mount Holder
DEHWM	Double Scope Wall Mount Holder

www.zutron.com



877-343-5873