

Zutron Medical would like to make your day easier with our **Multi-Enzymatic Cleansers**



Product Features

- **Highly effective formula**

These multi-enzymatic cleansers effectively clean away blood, proteins, fat, carbohydrates and synthetic liquids from surgical instruments.

- **Safe and easy to use**

Neutral pH formula is safe for all instruments, scopes and other healthcare equipment and is 100% biodegradable.

- **Cleansers are ideal for both manual and automated processes**

As well as ultrasonic units and washer sterilizers and decontaminators.

- **Non-Foaming formula**

Perfect for manual cleaning that demands high visibility and automated processes that require low foaming cleansers.

Multi-Enzymatic Cleanser Ordering Information

Order #	Description
ZUTR12071	Multi-Enzymatic Cleanser
ZUTR28601	Non-Foaming Multi-Enzymatic Cleanser
ZUTR50505	Concentrated Enzymatic Z-Sponge

Accessories Ordering Information

Order #	Description
ZUTR10ZDP	1 Ounce Dispensing Pump
ZUTRDEMA	Automatic Detergent Dosing Proportioner

Product Specifications

- Enzymatic Cleansers come in a case of 4, 1 gallon bottles.
 - *Case includes 1 ounce dispensing pump
 - *Diluted at 1/2 ounce to 1 ounce per gallon of water.
- Sponges come in a case of 80 sponges
 - *Diluted in up to 10 gallons of water

www.zutron.com



877-343-5873