


Zutron Medical is proud to announce the next generation **RTU Foamer 3.5**



Product Specifications

Dimensions:

15"(L) x 15"(W) x 31"(H)

Weight:

48 pounds

Style

Cart w/ 5 gallon pail holder

Electrical:

AC 110-120V, 2.5A, 60Hz

Warranty:

One year

Product Features

• Thick foaming power

Enzymatic foam prevents bioburden from drying on surgical instruments and keeps instruments moist for long periods of time, freeing up work flow.

• Easy to use

Fast application and powerful enzymes greatly reduce and ease the manual cleaning process, which results in quicker turn around time.

• New RTU Foamer 3.5 features

The new Fluid Control Knob allows user to adjust foam quality and greatly reduce the amount of chemical usage. Turn up for a more liquefied foamer or down for a airier foam that uses up to 45% less cleanser.

• Performance upgrades

Performance enhancements prolong life of the unit such as a cooling fan that helps air flow through the unit.

RTU Foamer 3.5 Ordering Information

Order #	Description
ZUTR50134	RTU Foamer Unit 3.5

RTU Foamer 3.5 Usage & Lease Program

Order #	Description
ZUTR50502-1	Free RTU Foamer 3.5 and 1 pail per month
ZUTR50502-2	Free RTU Foamer 3.5 and 2 pails per month
ZUTR50502-3	Free RTU Foamer 3.5 and 3 pails per month
ZUTR50502-4	Free RTU Foamer 3.5 and 4 pails per month

Accessory Ordering Information

Order #	Description
ZUTR50502	Multi-Enzymatic Foamer Cleanser RTU - 5 gal. pail
ZUTR50105	RTU Foamer 3.5 Replacement Wand
ZUTR50308	RTU Foamer 3.5 Replacement Nozzle
ZUTR50307	5 Gallon Pail Wrench
ZUTRMAT153	Super Absorbent Spill Mat Kit

www.zutron.com



877-343-5873